



Forty Years of
Swiss Innovation



LipoSurg

Liposuction PAL System with Integrated Infiltration Pump for TLA



Liposuction has never been easier!

It's time for the world innovative of the all new LipoSurg, a unit developed exactly to the needs of the modern Liposuction. Nouvag LipoSurg is the first of it's kind Liposuction unit combining two task in once - a high Flow TLA Infiltration Pump and a powerful reciprocation Liposuction Cannula.






Forty Years of
Swiss Innovation



Advantages / Reason Why

- ✓ Best price / performance ratio
- ✓ integrated Infiltration Pump
- ✓ High performance Power Cannula with stroke frequency from 3'300 to 5000/min
- ✓ clearly arranged display with all information at a glance
- ✓ Powerful Infiltration pump up to 450ml/min – tissue protecting
- ✓ lightweight and ergonomically handle for Power Cannula
- ✓ Safety first – Integrated Infiltration Pump to prevent damages or manipulation on pump

System overview



Power Cannula Handle

- ✓ Fast and easy suction of the fat with fast forward/backward movements up to 5'000 strokes /min
- ✓ Low vibration and lightweight handles for tiring free operating
- ✓ ergonomically Handle with High grade material
- ✓ enlarged suction channel for optimal extraction of the fat



Perfectly to combine with



Typ:	Vacuson 60
Betriebsspannung:	115/230 Volt bei 50-60 Hz
I-max.:	250 W
Saugleistung:	Bis 60 l/min.
Max. Vakuum:	-0,9 bar bei 686 mmHG
Abmessungen (WxDxH):	360 x 290 x 280
Gewicht:	12 kg

Antistatic wheels -2 breakable 3 Power outlets

Infiltration accessories

Cannulae for Infiltration

Ø Cannula	Length	Art.-No.
3.0 mm	250 mm	4350



Needle for Infiltration

Ø Needle	Length	Art.-No.
0.9 mm	70 mm	4561
1.0 mm	120 mm	4562



Packing unit with 100 sterile pieces.

Adapter with Luer-Lock attachment, to use with Cannulae with handle connector to Luer-Lock Injection needles. **Art. No.: 4398**



Infiltration of Tumescant Local Anaesthetic Solution



Luer-Lock Tube connection to Infiltration Pump.



Liposuction cannulas



Straight Cannulae for Abdominal Surgery

Ø Cannula	Length	Holes	Ø of Holes	Art.-No.
1.5 mm	150 mm	1	oval	4361
2.0 mm	150 mm	1	oval	4364
2.0 mm	150 mm	18	1.0 mm	4373
3.0 mm	150 mm	18	1.5 mm	4374
3.0 mm	200 mm	22	1.5 mm	4378
3.0 mm	300 mm	30	1.5 mm	4387
4.0 mm	200 mm	22	2.0 mm	4379
4.0 mm	300 mm	30	1.5 mm	4388

Angled Cannulae for Thigh Liposuction

Ø Cannula	Length	Holes	Ø of Holes	Art.-No.
3.0 mm	200 mm	22	1.5 mm	4381

Bowed Cannulae for Thigh Liposuction

Ø Cannula	Length	Holes	Ø of Holes	Art.-No.
3.0 mm	200 mm	22	1.5 mm	4362
3.0 mm	300 mm	30	1.5 mm	4365
4.0 mm	200 mm	22	1.5 mm	4368
4.0 mm	300 mm	30	1.5 mm	4372



Ø Cannula	Length	Art.-No.
3.0 mm	200 mm	4480
3.0 mm	300 mm	4481
4.0 mm	300 mm	4482
5.0 mm	300 mm	4483
6.0 mm	300 mm	4484



Ø Cannula	Length	Art.-No.
2.0 mm	150 mm	4364
3.0 mm	200 mm	4485
3.0 mm	300 mm	4486



Mercedes

Ø Cannula	Length	Art.-No.
3.0 mm	150 mm	4487
3.0 mm	250 mm	4488
4.0 mm	150 mm	4489
4.0 mm	250 mm	4490




Toledo

Ø Cannula	Length	Art.-No.
2.0 mm	50 mm	4595
3.0 mm	150 mm	4596
4.0 mm	250 mm	4597
5.0 mm	300 mm	4598



NouClean Maintenance Spray

- ✓ For maintenance of your top-quality NOUVAG precision instruments
- ✓ Environmentally friendly - without CFCs
- ✓ Reusable nozzles; easy change-over to full spray can
- ✓ Nozzle for contra angle
- ✓ Nozzle for micro motor 21ESS





Forty Years of
Swiss Innovation



General information

Nouvag guarantee:

- ✓ Maintain all relevant standards and directives
- ✓ Made out of high quality Medical stainless steel
- ✓ fully autoclave able and washable

