


Water baths WB / WDB / PB / HDS



Medical device class I


Heuser Apparatebau GmbH

Am Schlagbaum 10
D-42781 Haan

Postfach 3113
D-42769 Haan

Telefon +49 (0) 21 29 - 94 03 - 0
Telefax +49 (0) 21 29 - 94 03 - 33

www.heuser-haan.de
info@heuser-haan.de



Appropriation

Type WB: Warming and tempering of heat packs and bandages in physiotherapy

Type HDS: Damping of haypacks

Type WDB: Damping of thermoplastic bandages

Type PB: Warming and tempering of paraffin or fango-paraffin

Equipment

- Housing made of stainless steel with structured surface
- Top lid with handle bar in red, blue, green, yellow or black
- Gas spring (except WB 4-30)
- Energy saving insulation
- Electronic temperature control with digital temperature display (except WB 4-30)
- WB 4-30 with thermostatic temperature control, no temperature display
- Temperature measurement in the water
- Silicone heating element with large surface for a fast and consistent heating process
- Power On indicator lamp
- Heating indicator lamp
- Water outlet on the front side, ball valve 1/2"
- Stainless steel insert grid

Optional Accessory

Mobile steel rack, consisting of

- Powder coated steel frame
- Four castor rolls (two of which are lockable)
- Four support points for the water bath
- Stainless steel storage shelf

Device	Number of heat packs	Heating capacity	Outside dimensions mm L x W x H	Inside dimensions mm L x W x H	Equipment
WB 4-30/4	5	900 W	575x365x370	500x300x200	Insert grid, ball valve
WB 5-30/4	5	900 W	575x365x355	500x300x200	Insert grid, ball valve, gas spring, insulated top lid
WB 6-50/4	8	1500 W	675x470x385	590x400x200	Insert grid with 3 register bars, ball valve, gas spring, insulated top lid
WB 8-90/4	12	1500 W	725x470x490	645x400x320	Insert grid with 4 register bars, ball valve, gas spring, insulated top lid
WB 16-130/4	18	1500 W	775x490x515	695x420x340	Insert grid with 7 register bars, ball valve, gas spring, insulated top lid
WDB 1/4	Thermoplastics	1500 W	770x430x370	685x360x185	Insert grid, ball valve, gas spring, insulated top lid
WDB 2/4	Thermoplastics	900 W	570x420x280	485x350x130	Insert grid, ball valve, gas spring, insulated top lid
WDB 5/4	Thermoplastics	3000 W	775x620x375	695x550x200	Insert grid, ball valve, gas spring, insulated top lid
WB 6-100/4	Thermoplastics	1500 W	675x470x250	590x400x100	Insert grid, ball valve, gas spring, insulated top lid
HDS 5-30/4	Haypacks	900 W	575x365x355	500x300x200	Insert grid, ball valve, gas spring, insulated top lid
HDS 6-50/4	Haypacks	1500 W	675x470x385	600x400x200	Insert grid with 3 register bars, ball valve, gas spring, insulated top lid
HDS 8-90/4	Haypacks	1500 W	725x470x490	645x400x320	Insert grid with 4 register bars, ball valve, gas spring, insulated top lid
HDS 16-130/4	Haypacks	1500 W	775x490x515	695x420x340	Insert grid with 7 register bars, ball valve, gas spring, insulated top lid
PB 5-30/4	Paraffin, Fango-paraffin	900 W	575x365x355	500x300x200	Paraffin container 10 l, ball valve, gas spring, insulated top lid
PB 6-50/4	Paraffin, Fango-paraffin	1500 W	675x470x385	590x400x200	Paraffin container 20 l, ball valve, gas spring, insulated top lid

Technical details and equipment are subject to alteration!